



DZHK

DEUTSCHES ZENTRUM FÜR
HERZ-KREISLAUF-FORSCHUNG E.V.

Alignment of Multi-Sensor Data: Adjustment of Sampling Frequency and Time Shifts

Marcus Vollmer

Institute of Bioinformatics, University Medicine Greifswald

Young-DZHK Retreat 2018

08 June 2018, Majuwie Greifswald

Outline

- 1 Multi-Sensor Data
- 2 Adjustment of Time Shifts
- 3 Correction of Sampling Frequency
- 4 Perspectives

1. Multi-Sensor Data

What is multi-sensor data?

Multisensor data fusion, PK Varshney

Electronics & Communication Engineering Journal, 1997

Multisensor data fusion refers to the acquisition, processing and synergistic combination of information gathered by various knowledge sources and sensors to provide a better understanding of a phenomenon.

- ▶ Improved system reliability and robustness
- ▶ Extended coverage
- ▶ Increased confidence
- ▶ Shorter response time
- ▶ Improved resolution


Research team




Research team



Experimental setup



Experimental setup




- ▶ 5 min standing rest
- ▶ 5 min walking on treadmill (1.2 m/s)
- ▶ cognitive test (2-back audio test)
- ▶ 5 min walking on treadmill (1.2 m/s, 15% gradient)

In between: NASA Task Load Index to measure individual strain


Experimental setup

- ▶ EMotion Faros 360° | 1000 Hz
- ▶ SOMNOtouch NIBP | 512 Hz
- ▶ NeXus-10 MKII | 8192 Hz
- ▶ Polar RS800 Multi | 1000 Hz
- ▶ Hexoskin Smart Shirt | 256 Hz




Experimental setup

- ▶ EMotion Faros 360° | 1000 Hz
- ▶ SOMNOtouch NIBP | 512 Hz
- ▶ NeXus-10 MKII | 8192 Hz 8000 Hz
- ▶ Polar RS800 Multi | 1000 Hz
- ▶ Hexoskin Smart Shirt | 256 Hz




2. Adjustment of Time Shifts


Different starting time



Different starting time




Different starting time




Different starting time



Different starting time



Adjustment using ECG annotations



3. Correction of Sampling Frequency

Why should I adjust sampling frequencies?

- ▶ User manual says 1000 Hz
Expect an imprecise clock signal (quartz crystal): **+0.01%**
In 24h-measurement we assume to collect $1000 \cdot 24 \cdot 60 \cdot 60 = 86.400.000$ samples
Actually given the impreciseness we observe **86.408.640** samples

 $8.640 / 1000 \text{ Hz} = 8.64 \text{ s}$

Why should I adjust sampling frequencies?

- ▶ User manual says 1000 Hz
Expect an imprecise clock signal (quartz crystal): **+0.01%**
In 24h-measurement we assume to collect $1000 \cdot 24 \cdot 60 \cdot 60 = 86.400.000$ samples
Actually given the impreciseness we observe **86.408.640** samples


 $8.640 / 1000 \text{ Hz} = 8.64 \text{ s}$
- ▶ Identical signals are not aligned anymore
- ▶ Without adjustment wrong conclusions can be concluded:

Why should I adjust sampling frequencies?

- ▶ User manual says 1000 Hz
Expect an imprecise clock signal (quartz crystal): **+0.01%**
In 24h-measurement we assume to collect $1000 \cdot 24 \cdot 60 \cdot 60 = 86.400.000$ samples
Actually given the impreciseness we observe **86.408.640** samples

 $8.640 / 1000 \text{ Hz} = 8.64 \text{ s}$
- ▶ Identical signals are not aligned anymore
- ▶ Without adjustment wrong conclusions can be concluded:
"The heart rate increased subsequent to the accident."
"The heart rate increased just before the accident."


Linear adjustment of sampling frequencies



- ▶ Pairwise differences s_i are linear decreasing or increasing with time





Linear adjustment of sampling frequencies

$$Fs_correct = Fs * (1 - \text{'Slope of robust regression fit'})$$



Is a linear adjustment sufficient?

Hexoskin as reference:



4. Perspectives

Benefits from multi-sensor data

Well aligned signals



Benefits from multi-sensor data

Well aligned signals



- ▶ Quality assurance
- ▶ Increased confidence
- ▶ Backup signals
(low signal quality, signal loss)
- ▶ Use of full signal variety
(body movement, temperature, respiration)

Benefits from multi-sensor data

Well aligned signals




- ▶ Quality assurance
- ▶ Increased confidence
- ▶ Backup signals
(low signal quality, signal loss)
- ▶ Use of full signal variety
(body movement, temperature, respiration)
- ▶ Comparison of manufacturers heart rate estimation
- ▶ Extraction of realistic noise
- ▶ Verification/accuracy of Polar RR Intervals

Benefits from multi-sensor data



► Verification/accuracy of Polar RR Intervals

Future experiments





Thank You for
Your Attention!